



Recombinant *Oryza sativa* 50S ribosomal protein L2, chloroplastic (rpl2-A)

Product Code	CSB-YP315020OCO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C495
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAKHLYKTPI PSTRKGTIDR QVKSNPRNNL IHGRHRCGKG RNSRGIITAR HRGGGHKRLY RKIDFRRNQK DISGRIVTIE YDPNRNAYIC LIHYGDGEKG YILHPRGAIL GDTIVSGTKV PISMGNALPL TDMP LGTAIH NIEITRGRGG QLARAAGAVA KLIAKEGKSA TLRLPSGEVR LVSQNCLATV GQVGNVGVNQ KSLGRAGSKC WLGKRPVVRG VVMNPVDHPH GGEGKAPIG RKKPTTPWGY PALGRRTRKR KKYSDFSILR RRK
Source	Yeast
Target Names	rpl2-A
Protein Names	Recommended name: 50S ribosomal protein L2, chloroplastic
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.