



# Recombinant Avian leukosis virus RSA Gag-Pro polyprotein (gag-pro), partial

<b>Product Code</b>	CSB-BP315187ASN
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P0C776
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Avian leukosis virus RSA (RSV-SRA) (Rous sarcoma virus (strain Schmidt-Ruppin A))
<b>Purity</b>	>85% (SDS-PAGE)
<b>Source</b>	Baculovirus
<b>Target Names</b>	gag
<b>Protein Names</b>	Recommended name: Gag-Pro polyprotein Cleaved into the following 8 chains: 1. Matrix protein p19 2. p2A 3. p2B 4. p10 5. p3 6. Capsid protein p27 7. Nucleocapsid protein p12 8. Protease p15 EC= 9. 3.4.23.-
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.