



Recombinant Avian leukosis virus RSA Gag-Pro polyprotein (gag-pro), partial

| | |
|--------------------------|---|
| Product Code | CSB-YP315187ASN |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P0C776 |
| Product Type | Recombinant Protein |
| Immunogen Species | Avian leukosis virus RSA (RSV-SRA) (Rous sarcoma virus (strain Schmidt-Ruppin A)) |
| Purity | >85% (SDS-PAGE) |
| Source | Yeast |
| Target Names | gag |
| Protein Names | Recommended name: Gag-Pro polyprotein Cleaved into the following 8 chains: 1. Matrix protein p19 2. p2A 3. p2B 4. p10 5. p3 6. Capsid protein p27 7. Nucleocapsid protein p12 8. Protease p15 EC= 9. 3.4.23.- |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |