



Recombinant Escherichia coli 3-oxoadipyl-CoA/3-oxo-5,6-dehydrosuberil-CoA thiolase (paaJ)

Product Code	CSB-EP315200ENV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C7L2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MREAFICDGI RTPIGRYGGA LSSVRADDLA AIPLRELLVR NPRLDAECID DVILGCANQA GEDNRNVARM ATLLAGLPQS VSGTTINRLC GSGLDALGFA ARAIKAGDGD LLIAGGVESM SRAPFVMGKA ASAFSRQAEM FDTTIGWRFV NPLMAQQFGT DSMPETAENV AELLKISRED QDSFALRSQQ RTAKAQSSGI LAEIIVPVVL KNKKGVVTEI QHDEHLRPET TLEQLRGLKA PFRANGVITA GNASGVNDGA AALIIASEQM AAAQGLTPRA RIVAMATAGV EPRLMGLGPV PATRRVLERA GLSIHMDMDVI ELNEAFAAQA LGVLRGLLP DDAPHVNPNG GAIALGHPLG MSGARLALAA SHELHRRNGR YALCTMCIGV GGGIAMILER V
Source	E.coli
Target Names	paaJ
Protein Names	Recommended name: 3-oxoadipyl-CoA/3-oxo-5,6-dehydrosuberil-CoA thiolase EC= 2.3.1.174 EC= 2.3.1.n11
Expression Region	1-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.