



# Recombinant Human Putative protein FAM200B (FAM200B)

<b>Product Code</b>	CSB-EP315397HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P0CF97
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDHFFIKRKR NSEVKYTEAC SSSSVESGIV NSDNIEKNTD SNLQTSTSFE PHFKKKKVSA RRYNEDYLKY GFIKCEKPF E NDRPQCVICN NILANESLKP SKLKRHLETQ HAELIDKPLE YFQRKKKDIK LSTQFLSCST AVSEKALLSS YLVAYRVAKE KiantAAEKI ILPACLDMVR TIFDDKSADK LKTIPNDNTV SLRICTIAEH LETMLITRLQ SGIDFAIQLD ESTDIGSCTT LLVYVRYAWQ DDFLEDFLCF LNLTSHLSGL DIFTELERRI VGQYKLNWKN CKGITS SDGTA TMTGKHSRVI KKLLEVTNNG AVWNHCFIHR EGLASREIPQ NLMEVLKNAV KVVNFIKSS LNSRLLETFC SEIGTNHHL LYHTKIRWLS QGKILSRVYE LRNEIHFFLI EKKSHLASIF EDDTWVTKLA YLTDIFSILN ELSLKLQGKN SDVFQHVERI QGFRKTL LLW QVRLKSNRPS YYMFPRFLQH IEENIINENI LKEIKLEILL HLTSLSQTFN HFFPEEKFET LRENSWVKDP FAFRHPESII ELNLVPEEEN ELLQLSSSYT LKNDYETLSL SAFWMKV KED FPLLSRKS VL LLL PFTTSL CELGFSILTQ LKTKERNGLN CAAVMRVALS SCVPDWNELM NRQAHPS
<b>Source</b>	E.coli
<b>Target Names</b>	FAM200B
<b>Protein Names</b>	Recommended name: Putative protein FAM200B
<b>Expression Region</b>	1-657
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.