



Recombinant Human Putative protein FAM200B (FAM200B)

Product Code	CSB-MP315397HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CF97
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDHFFIKRKR NSEVKYTEAC SSSSVESGIV NSDNIEKNTD SNLQTSTSFE PHFKKKKQVSA RRYNEDYLKY GFIKCEKPF E NDRPQCVCIN NILANESLKP SKLKRHLETQ HAELIDKPLE YFQRKKKDIK LSTQFLSCST AVSEKALLSS YLVAYRVAKE KiantAAEKI ILPACLDMVR TIFDDKSADK LKTIPNDNTV SLRICTIAEH LETMLITRLQ SGIDFAIQLD ESTDIGSCTT LLVYVRYAWQ DDFLEDFLCF LNLTSHLSGL DIFTELERRI VGQYKLNWKN CKGITS DGTA TMTGKHSRVI KKLLEVTNNG AVWNHCFIHR EGLASREIPQ NLMEVLKNAV KVVNFIKSS LNSRLLETFC SEIGTNHHL LYHTKIRWLS QGKILSRVYE LRNEIHFFLI EKKSHLASIF EDDTWVTKLA YLTDIFSILN ELSLKLQGN SDVFQHVERI QGFRKTL LLW QVRLKSNRPS YYMFP RFLQH IEENIINENI LKEIKLEILL HLTSLSQTFN HFFPEEKFET LRENSWVKDP FAFRHPESII ELNLVPEEEN ELLQLSSSYT LKNDYETLSL SAFWMKV KED FPLLSRKS VL LLLPTTTSL CELGFSILTQ LKTKERNGLN CAAVMRVALS SCVPDWNELM NRQAHPS
Source	Mammalian cell
Target Names	FAM200B
Protein Names	Recommended name: Putative protein FAM200B
Expression Region	1-657
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.