



Recombinant D-tagatose-1,6-bisphosphate aldolase subunit gatZ (gatZ)

Product Code	CSB-EP315419ENL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C8J9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MKTLIARHKA GEHIGICSVC SAHPLVIEAA LAFDRNSTRK VLIETSQV NQFGGYTGMT PADREFVFT IADKVGFARE RIILGGDHLG PNCWQQENAD AAMEKSVELV KAYVRAGFSK IHLDASMSCA GDPIPLAPET VAERAAVLCF AAESVATDCQ REQLSYVIGT EVPVPGGEAS AIQSVHITRV EDAANTLRTH QKAFIARGLA EALTRVIAIV VQPGVEFDHS NIIHYQPQEA QPLAQWIENT RMVYEAHSTD YQTRTAYWEL VRDHFAILKV GPALTFALRE AIFALAQIEQ ELIAPENRSG CLAVIEEVMF DEPQYWKKYY RTGFNDLLD IRYSLSDRIR YYWPHSRIKN SVETMMVNLE GMEIPLGMIS QYLPKQFERI QSGELSAIPH QLIMDKIYDV LRAYRYGCAE
Source	E.coli
Target Names	gatZ
Protein Names	Recommended name: D-tagatose-1,6-bisphosphate aldolase subunit gatZ
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.