



Recombinant Zea mays Stress-related protein 1 (SRP1)

Product Code	CSB-YP315420ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C8M9
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MTSESSTPTG STRALPASIT RSSSSTLSTR PSASTPASTG SLTPSTSTAS LVVSPPPARS PVVSTRATAT TRPRLAAERP GSATTPCPCG DTDKQLAPCL QEGAGFVIRN AAVRLYMRFG RVVFWHCMLG LWGVGNLIW
Source	Yeast
Target Names	SRP1
Protein Names	Recommended name: Stress-related protein 1
Expression Region	1-139
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.