



# Recombinant Human Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial (PDHA1), partial

<b>Product Code</b>	CSB-MP017715HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P08559
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ANDATFEIKK CDLHRLEEGP PVTTVLTRD GLKYYRMMQT VRRMELKADQ LYKQKIIRGF CHLCDGQEAC CVGLEAGINP TDHLITAYRA HGFTFTRGLS VREILAEITG RKGCGAKGKG GSMHMYAKNF YGGNGIVGAQ VPLGAGIALA CKYNGKDEVC LTLYGDGAAN QGQIFEAYNM AALWKLPCIF ICENNRYGMG TSVERAAAST DYYKRGDFIP GLRVDGMDIL CVREATRFAA AYCRSGKGPI LMELQTYRYH GHSMSDPGVS YRTREEIQEV RSKSDPIMLL KDRMVNSNLA SVEELKEIDV EVRKEIEDAA QFATADPEPP LEELGYHIYS SDPPFEVRGA NQWIKFKSVS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PDHA1
<b>Protein Names</b>	Recommended name: Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial EC= 1.2.4.1Alternative name(s): PDHE1-A type I
<b>Expression Region</b>	31-390
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.