



# Recombinant Escherichia coli Transposase insH for insertion sequence element IS5H (insH6)

<b>Product Code</b>	CSB-MP315542ENV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P0CE53
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFVIWSHG TG FIMSHQLTFA DSEFSSKRRQ TRKEIFLSRM EQILPWQNMV EVIEPFYPKA GNGRRPYPLE TMLRIHCMQH WYNLSDGAME DALYEIASMR LFARLSLDSA LPDRTTIMNF RHLLEQHQLA RQLFKTINRW LAEAGVMMTQ GTLVDATIE APSSTKNKEQ QRDPEMHQTK KGNQWHFGMK AHIGVDAKSG LTHSLVTTAA NEHDLNQLGN LLHGEEQFVS ADAGYQGAPQ REELAEVDVD WLIAERPGKV RTLKQHPRKN KTAINIEYMK ASIRAKVEHP FRIIKRQFGF VKARYKGLLK NDNQLAMLFT LANLFRADQM IRQWERSH
<b>Source</b>	Mammalian cell
<b>Target Names</b>	insH6
<b>Protein Names</b>	Recommended name: Transposase insH for insertion sequence element IS5H
<b>Expression Region</b>	1-338
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.