



Recombinant Escherichia coli Transposase insH for insertion sequence element IS5-18 (JW5951)

Product Code	CSB-MP315544ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CE65
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MSHQLTFADS EFSSKRRQTR KEIFLSRMEQ ILPWQNMVEV IEPFYPKAGN GRRPYPLETM LRIHCMQHWY NLSDGAMEDA LYEIASMRLF ARLSLDSALP DRTTIMNFRH LLEQHQLARQ LFKTINRWLA EAGVMMTQGT LVDATIIEAP SSTKNKEQQR DPEMHQTKKG NQWHFGMKAH IGVDKASGLT HSLVTTAANE HDLNQLGNLL HGEEQFVSAD AGYQGAPQRE ELAEVDVDWL IAERPGKVRT LKQHPRKNKT AINIEYMKAS IRARVEHPFR IIKRQFGFVK ARYKGLLKND NQLAMLFTLA NLFRADQMIR QWERSH
Source	Mammalian cell
Target Names	JW5951
Protein Names	Recommended name: Transposase insH for insertion sequence element IS5-18
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.