



Recombinant Escherichia coli Transposase insF for insertion sequence IS3B (insF2)

Product Code	CSB-BP315560ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CF80
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKYVFIEKHQ AEFSIKAMCR VLRVARSGWY TWCQRRTRIS TRQQFRQHCD SVVLAAFTRS KQRYGAPRLT DELRAQGYPF NVKTVAASLR RQGLRAKASR KFSPVSYRAH GLPVSENLE QDFYASGPNQ KWAGDITYLR TDEGWLYLAV VIDLWSRAVI GWSMSPRMTA QLACDALQMA LWRRKRPRNV IVHTDRGGQY CSADYQAQLK RHNLRGMSMA KGCCYDNACV ESFFHSLKVE CIHGEHFISR EIMRATVFN Y IECDYNRWRR HSWCGGLSPE QFENKNLA
Source	Baculovirus
Target Names	insF2
Protein Names	Recommended name: Transposase insF for insertion sequence IS3B
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.