



Recombinant Escherichia coli Putative transposase InsL for insertion sequence element IS186B (insL2)

| | |
|--------------------------|--|
| Product Code | CSB-MP315562ENV |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P0CF92 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli (strain K12) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MNYSHDNWSA ILAHIGKPEE LDTSARNAGA LTRRREIRDA ATLLRLGLAY GPGGMSLREV TAWAQLHDVA TLSDVALLKR LRNAADWFGI LAAQTLAVRA AVTGCTSGKR LRLVDGTAIS APGGGSAEWR LHMGYDPHTC QFTDFELTDS RDAERLDRFA QTADEIRIAD RGFGSRPECI RSLAFGEADY IVRVHWRGLR WLTAEGMRFD MMGFLRGLDC GKNGETTVMI GNSGNKKAGA PFPARLIAVS LPPEKALISK TRLLSENRRK GRVVQAETLE AAGHVLLLS LPEDEYSAEQ VADCYRLRWQ IELAFKRLKS LLHLDALRAK EPELAKAWIF ANLLAAFLID DIIQPSLDFP PRSAGSEKKN |
| Source | Mammalian cell |
| Target Names | insL2 |
| Protein Names | Recommended name: Putative transposase InsL for insertion sequence element IS186B |
| Expression Region | 1-370 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |