



Recombinant *Saccharomyces cerevisiae* DNA damage-inducible protein 2 (DDI2)

Product Code	CSB-BP315583SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CH63
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSQYGFVRVP REVEKAIPVV NAPRPRAVVP PPNSETARLV REYAAKELTA PVLNHSRLRVF QYSVAIIRDQ FPAWDLQEV LYVTCLLDHI ATTDKNMRAT KMSFEYYGGI LSRELVFNAT GGNQDYADAV TEAIRHQDL TGTGYITTLG LILQIATTLN NVGSNTDLIH IDTVSAINEQ FPRLHWLSCF ATVVDTENS KRWGHTSSLG DDFSCKVICN TFGYN
Source	Baculovirus
Target Names	DDI2
Protein Names	Recommended name: DNA damage-inducible protein 2
Expression Region	1-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.