



Recombinant *Cryptococcus neoformans* var. *neoformans* serotype D Defective in cullin neddylation protein 1 (DCN1)

Product Code	CSB-MP315677CTK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CN07
Product Type	Recombinant Protein
Immunogen Species	<i>Cryptococcus neoformans</i> var. <i>neoformans</i> serotype D (strain B-3501A) (<i>Filobasidiella neoformans</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MASIVSPLST KDSLLVTQFR AITGTSSAEA AKYIKRYKHI EAAVDAFYNN EPAPRADPAQ ERKLGEIWEK FKDPSPDKLI KIDGTMELCE ELDIDPGTDA VLFCLAADLG SKATGEWEKA PFVAGIASYP GNIDSLPKLK AYLPTLRKKL VSDPEYFKKV YNHAFQLARG GPQSLTRSLP LDTAIDLWTL FFPPAFNHSP SALSHLPDMS PPQFTQPEFD LWIEFMQQKN KAVSKDTWAL LVDFARGIDK DFKEYDEDGA WPSMIDDFVE YVREQKRGQ
Source	Mammalian cell
Target Names	DCN1
Protein Names	Recommended name: Defective in cullin neddylation protein 1
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.