



Recombinant Human Testis-specific Y-encoded protein 8 (TSPY8)

Product Code	CSB-YP315760HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CW00
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRPEGSLTYW VPERLRQGFC GVGRAAQALV CASAKEGTAF RMEAVQEGAA GVESEQAALG EEAVLLDDI MAEEVVAAE EGLVERREEA QRAQAVPGP GPMTPEALE ELLAVQVELE PVNAQARKAF SRQREKMERR RKPHLDRRGA VIQSVPGFWA NVIANHPQMS ALITDEDEDM LSYMVSLEVE EEKHPVHLCK IMLFFRSNPY FQNKVITKEY LVNITEYRAS HSTPIEWYLD YEVEAYRRRH HNSSLNFFNW FSDHNFAGSN KIAEILCKDL WRNPLQYYKR MKPPEEGTET SGDSQLLS
Source	Yeast
Target Names	TSPY8
Protein Names	Recommended name: Testis-specific Y-encoded protein 8
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.