



# Recombinant *Borrelia burgdorferi* Basic membrane protein D (bmpD)

<b>Product Code</b>	CSB-YP316105BUD
<b>Storage</b>	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
<b>Uniprot No.</b>	P0CL55
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Borrelia burgdorferi</i> (strain ATCC 35210 / B31 / CIP 102532 / DSM 4680)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	CSSS DDGKSEAKTV SLIVDGAFDD KGFNESSSKA IRCLKADLNI NIEKASTGN SYLGDIANLE DGNSNLIWGI GFRLSDILFQ RASENVSVNY AIIIEGVYDEI QIPKNLLNIS FRSEEVAFLA GYFASKASKT GKIGFVGGVR GKVLESFMYG YEAGAKYANS NIKVVSQYVG TFGDFGLGRS TASNMYRDGV DIIFAAAGLS GIGVIEAAKE LGPDHYIIGV DQDQSYLAPN NVIVSAVKKV DSLMYSLTKK YLETGVLDGG KTMFLGLKED GLGLVLNENL KSNYSEIYNK SLKIGQSIMN GIIKVPYDKV SYDNFVLQME N
<b>Source</b>	Yeast
<b>Target Names</b>	bmpD
<b>Protein Names</b>	Recommended name: Basic membrane protein D
<b>Expression Region</b>	17-341
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.