



Recombinant *Borrelia burgdorferi* Basic membrane protein C (bmpC)

Product Code	CSB-MP316108BUD
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CL65
Product Type	Recombinant Protein
Immunogen Species	<i>Borrelia burgdorferi</i> (strain ATCC 35210 / B31 / CIP 102532 / DSM 4680)
Purity	>85% (SDS-PAGE)
Sequence	CFKS NKKSIKSDKV VGVLAHGSF YDKGYNQSVH DGVVKLRDNF GIKLITKSLR PYPPIEGKRL TVDEAMTEDA YEVQKNPLNL FWLIGYRFSD LSVKLSYERP DIYYGIIDAF DYGDIQVPKN SLAIKFRNEE AAFLAGYIAA KMSRKEKIGF LTGPMSEHVK DFKFGFKAGI FYANPKLRLV SKKAPSLFDK EKGKAMALFM YKEDKVGIVF PIAGITGLGV YDAAKELGPK YYVIGLNQDQ SYIAPQNVIT SIIKDIGKVI YSISSEYINN RVFKGGIIID RGLKEGVIEI VKDPDVLNNR LVDEVIDLEN KIISGEIIVP DSEYAFDLFK SKL
Source	Mammalian cell
Target Names	bmpC
Protein Names	Recommended name: Basic membrane protein C
Expression Region	17-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.