



Recombinant *Aspergillus clavatus* Ribonuclease clavin (cla)

Product Code	CSB-BP316109AVC
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CL71
Product Type	Recombinant Protein
Immunogen Species	<i>Aspergillus clavatus</i> (strain ATCC 1007 / CBS 513.65 / DSM 816 / NCTC 3887 / NRRL 1)
Purity	≥85% (SDS-PAGE)
Sequence	AAT WTCMNEQKNP KTNKYENKRL LYNQNNAESN AHHAPLSDGK TGSSYPHWFT NGYDGDGKIL KGRTPIKWGN SDCDRPPKHS KNGDGKNDHY LLEFPTFPDG HQYNFDSKPP KEDPGPARVI YTYPNKVFCG IVAHTRENQG DLKLCSH
Source	Baculovirus
Target Names	cla
Protein Names	Recommended name: Ribonuclease clavin EC= 3.1.27.-
Expression Region	28-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.