



Recombinant *Cryptococcus neoformans* var. *neoformans* serotype D Peptidyl-prolyl isomerase CWC27 (CWC27)

Product Code	CSB-BP316158CTL
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CP92
Product Type	Recombinant Protein
Immunogen Species	<i>Cryptococcus neoformans</i> var. <i>neoformans</i> serotype D (strain JEC21 / ATCC MYA-565) (<i>Filobasidiella neoformans</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSNLYATEPA TNGKVIIDTT AGEIEVELWG KECPKAVRNF LALTMEGYYD GVIFHRVVPF FIIQSGDPTG TGMGGESFYG EPFEDEIHGR LKFNRRGLLG MANNGSRNSN TSQFFITLDA APELTNKHTM FGKIVGNTIF NVLNIGNLDT DKEERPLIPP KIRRIHIIEN PFDDIVPRIT AEEKKAQQKA KLEAKKDMEQ RERRAKAKKN TGLLSFGDSE EIPEEEVTVK KKSMTQRDLV DPSEAEHTKS KTSKMTETFV NIPPSLKDLG KSRENEASEE KKAVDLKNIR AQHEREKAGG SAARQAEIKR MEEDLRRLKK RSGSVSDSES DSSSRARRKG PSYLEQELAK YASKRGRAAM KAGNKRGRRD EEEDVLTEMR KFSKRVMQAG DEPEEEQAEE IEEGEAKEEG TGIGAAMAE EGGIEVDDDV GWLTHKLFQ VDDKELTRRA EDEYAVIDPR AKARDLLGKP DKKKLGKGNPN RRTVRNSGRN R
Source	Baculovirus
Target Names	CWC27
Protein Names	Recommended name: Peptidyl-prolyl isomerase CWC27 Short name= PPlase CWC27 EC= 5.2.1.8 Alternative name(s): Rotamase CWC27
Expression Region	1-491
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.