



Recombinant Human Coiled-coil domain-containing protein 166 (CCDC166)

Product Code	CSB-MP316213HU
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CW27
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAPKKKRGPS AGSQPGGAAA AGAEQPLSER AQYLQREHAL LSEQLDTCEE SVDQVLRENA FLDREALRLR EENRLYASYV SARAQRCACA IVRLDEQNRV DLAQIHWQRA ELASLYHGRE DGVRAQLLEM EARAAQMAQQ VQELQPYKVL QLEQLARIRA LERELLHMRV EHTQLLHRVK RRFLEDKAAF EREARQRVQS LARRAEREAV RALVAHTQAI KADNGRLRQE LLLLLRRTQL LHHTRRQLE QREQLHREHE DTRDLARVHG WLRRGPGGPP LWERPAFSQP TSRPGSLAAP ISPSRAASQT PSVVPSRAAP RASSVPSRE ASRVPSLVLS SMDSRVPSLA TSKVGSRMPS LTASRAGSRA LSLVQSLEGS GISSGSSPRV SSQDTLRSTK SGPKLLSGLS RDRDPALLPP QSEDSVNAEA AAEASPGRA
Source	Mammalian cell
Target Names	CCDC166
Protein Names	Recommended name: Coiled-coil domain-containing protein 166
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.