



Recombinant *Saccharomyces cerevisiae* 60S ribosomal protein L19-B (RPL19B)

Product Code	CSB-YP316244SVG
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CX83
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	ANLRTQKRL AASVVGVGKR KVLDPNETS EIAQANSRNA IRKLVKNGTI VKKAVTVHSK SRTRAHAQSK REGRHSGYGK RKGTREARLP SQVWVIRRLR VLRLLAKYR DAGKIDKHLY HVLYKESKGN AFKHKRALVE HIIQAKADAQ REKALNEEAE ARRLKNRAAR DRRAQRVAEK RDALLKEDA
Source	Yeast
Target Names	RPL19B
Protein Names	Recommended name: 60S ribosomal protein L19-B Alternative name(s): L23 RP15L RP33 YL14
Expression Region	2-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.