



Recombinant Mouse Transcription factor jun-B (Junb)

Product Code	CSB-EP011973MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P09450
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MCTKMEQPFY HDDSYAAAGY GRSPGSLSLH DYKLLKPTLA LNLADPYRGL KGPGARGPGP EGSGAGSYFS GQGS DTGASL KLASTELERL IVPNSNGVIT TTPTPPGQYF YPRGGGSGGG TGGGVTEEQE GFADGFVKAL DDLHKMNHVT PPNVSLGASG GPQAGPGGVY AGPEPPPVYT NLSSYSPASA PSGGSGTAVG TGSSYPTATI SYLPHAPPFA GGHPAQLGLS RGASAFKEEP QTVPEARSRD ATPPVSPINM EDQERIKVER KRLRNRLAAT KCRKRKLERI ARLEDKVKTL KAENAGLSSA AGLLREQVAQ LKQKVMTHVS NGCQLLLGVK GHAF
Source	E.coli
Target Names	Junb
Protein Names	Recommended name: Transcription factor jun-B Alternative name(s): MyD21
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.