



Recombinant *Saccharomyces cerevisiae* Putative formate dehydrogenase 2 (FDH2)

Product Code	CSB-MP316565SVG
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CF35
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVVINKQLMV SGILPAWLKN EYDLEDKIIS TVGAGRIGYR VLRLVAFNP KKLLYDYQE LPAEAINRLN EASKLFNGRG DIVQRVEKLE DMVAQSDVVT INCPLHKDSR GLFNKKLISH MKDGAYLVNT ARGaicVAED VAEAVKSGKL AGYGGDVWDK QPAPKDHPWR TMDNKDHVGN AMTVHISGTS LHAQKRYAQG VKNILNSYFS KKFdyRPQDI IVQNGSYATR AYGQKK
Source	Mammalian cell
Target Names	FDH2
Protein Names	Recommended name: Putative formate dehydrogenase 2 EC= 1.2.1.2 Alternative name(s): NAD-dependent formate dehydrogenase 2
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.