



Recombinant Escherichia coli Transposase insC for insertion element IS2-10 (JW5919)

Product Code	CSB-EP316567ENV
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CF49
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MIDVLGPEKR RRRTTQEKIA IVQQSFEPGM TVSLVARQHG VAASQLFLWR KQYQEGSLTA VAAGEQVPA SELAAAMKQI KELQRLLGKK TMENELLKEA VEYGRAKKWI AHAPLLPGDG E
Source	E.coli
Target Names	JW5919
Protein Names	Recommended name: Transposase insC for insertion element IS2-10
Expression Region	1-121
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.