



Recombinant Mouse Cytotoxic T-lymphocyte protein 4 (Ctla4), partial

Product Code	CSB-MP006163MO1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P09793
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	EAIQV TQPSVVLASS HGVASFPC EY SPSHNTDEV R VTVLRQTNDQ MTEVCATTF T EKNTVGFLDY PFCSGTFNES RVNLTIQGLR AVDTGLYLCK VELMYPPPYF VGMGNGTQIY VIDPEPCPDS D
Source	Mammalian cell
Target Names	Ctla4
Protein Names	Recommended name: Cytotoxic T-lymphocyte protein 4 Alternative name(s): Cytotoxic T-lymphocyte-associated antigen 4 Short name= CTLA-4 CD_antigen= CD152
Expression Region	36-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This gene is a member of the immunoglobulin superfamily and encodes a protein which transmits an inhibitory signal to T cells. The protein contains a V domain, a transmembrane domain, and a cytoplasmic tail. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. The membrane-bound isoform functions as a homodimer interconnected by a disulfide bond, while the soluble isoform functions as a monomer. Mutations in this gene have been associated with insulin-dependent diabetes mellitus, Graves disease, Hashimoto thyroiditis, celiac disease, systemic lupus erythematosus, thyroid-associated orbitopathy, and other autoimmune diseases.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.