



# Recombinant Escherichia coli Transposase insH for insertion sequence element IS5-15 (JW5928)

<b>Product Code</b>	CSB-EP316718ENV-B
<b>Storage</b>	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
<b>Uniprot No.</b>	P0CE62
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSHQLTFADS EFSSKRRQTR KEIFLSRMEQ ILPWQNMVEV IEPFYPKAGN GRRPYPLETM LRIHCMQHWY NLSDGAMEDA LYEIASMRLF ARLSLDSALP DRTTIMNFRH LLEQHQLARQ LFKTINRWLA EAGVMMTQGT LVDATIIEAP SSTKNKEQQR DPEMHQTKKG NQWHFGMKAH IGVDKSGLT HSLVTTAANE HDLNQLGNLL HGEEQFVSAD AGYQGAPQRE ELAEVDVDWL IAERPGKVRT LKQHPRKNKT AINIEYMKAS IRARVEHPFR IIKRQFGFVK ARYKGLLKND NQLAMLFTLA NLFRADQMIR QWERSH
<b>Source</b>	E.coli
<b>Target Names</b>	JW5928
<b>Protein Names</b>	Recommended name: Transposase insH for insertion sequence element IS5-15
<b>Expression Region</b>	1-326
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.