



Recombinant *Saccharomyces cerevisiae* 60S ribosomal protein L20-B (RPL20B)

Product Code	CSB-EP316941SVG-B
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CX24
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAHFKEYQVI GRRLPTESVP EPKLFMRMRF ASNEVIAKSR YWYFLQKLHK VKKASGEIVS INQINEAHPT KVKNFGVWVR YDSRSGTHNM YKEIRDVSRV AAVETLYQDM AARHRARFRS IHILKVAEIE KTADVQRQYV KQFLTKDLKF PLPHRVQKST KTF SYKRPST FY
Source	E.coli
Target Names	RPL20B
Protein Names	Recommended name: 60S ribosomal protein L20-B Alternative name(s): L18a
Expression Region	1-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.