



Recombinant Rhodobacter capsulatus Electron transport complex protein rnfB (rnfB)

Product Code	CSB-MP316977RLC
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CZ14
Product Type	Recombinant Protein
Immunogen Species	Rhodobacter capsulatus (Rhodopseudomonas capsulata)
Purity	≥85% (SDS-PAGE)
Sequence	MIAAAASMSA LGLGLGYLLG AAARKFHVET PPIVEEIAKI LPGTNCGACG FPGCNGLAEA MAEGNAPVTA CTPGGRDVAL ALAEIVTVEA GADAGPIAEI EPMVAFVFED HCTGCQKCFK RCPTDAIVGG AKQIHTVVM D ACIGCDACIE VCPT EAI VSR VKPKTLKTWY WDKPQPGLVA ASAETAA
Source	Mammalian cell
Target Names	rnfB
Protein Names	Recommended name: Electron transport complex protein rnfB Alternative name(s): Nitrogen fixation protein rnfB
Expression Region	1-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.