



Recombinant Horse Phospholipase A2 (PLA2G1B)

Product Code	CSB-MP018090HO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P00594
Product Type	Recombinant Protein
Immunogen Species	Equus caballus (Horse)
Purity	>85% (SDS-PAGE)
Sequence	AVW QFRSMIQCTI PNSKPYLEFN DYGCYCGLGG SGTPVDELDA CCQVHDNCYT QAKELSSCRF LVDNPYTESY KFSCSGTEVT CSDKNNACEA FICNCDRNAA ICFSKAPYNP ENKNLDSKRK CA
Source	Mammalian cell
Target Names	PLA2G1B
Protein Names	Recommended name: Phospholipase A2 EC= 3.1.1.4 Alternative name(s): Group IB phospholipase A2 Phosphatidylcholine 2-acylhydrolase 1B Cleaved into the following chain: 1. Phospholipase A2 isoform 2
Expression Region	8-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.