



Recombinant Pongo abelii NADH dehydrogenase [ubiquinone] iron-sulfur protein 8, mitochondrial (NDUFS8)

| | |
|--------------------------|---|
| Product Code | CSB-BP317146PYX |
| Storage | The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C. |
| Uniprot No. | P0CB97 |
| Product Type | Recombinant Protein |
| Immunogen Species | Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii) |
| Purity | >85% (SDS-PAGE) |
| Sequence | TYKYVN MQPEMDMKS VTDRAARTLL LTELFRGLGM TLSYLFREPA TINYPFKGP LSPFRGEHA LRRYPSGEER CIACKLCEAI CPAQAITIEA EPRADGSRRT TRYDIDMTKC IYCGFCQEAC PVDAIVEGPN FEFSTETHEE LLYNKEKLLN NGDKWEAEIA ANIQADYLYR |
| Source | Baculovirus |
| Target Names | NDUFS8 |
| Protein Names | Recommended name: NADH dehydrogenase [ubiquinone] iron-sulfur protein 8, mitochondrial EC= 1.6.5.3 EC= 1.6.99.3 Alternative name(s): Complex I-23kD Short name= CI-23kD NADH-ubiquinone oxidoreductase 23 kDa subunit |
| Expression Region | 35-210 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |