



Recombinant Escherichia coli Transposase insF for insertion sequence IS3fA (insF6)

Product Code	CSB-MP317178ENV
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CF84
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKYVFIEKHQ AEFSIKAMCR VLRVARSGWY TWCQRRTRIS TRQQFRQHCD SVVLAAFTRS KQRYGAPRLT DELRAQGYPF NVKTVAAASLR RQGLRAKASR KFSPVSYRAH GLPVSLENLLE QDFYASGPNQ KWAGDITYLR TDEGWLYLAV VIDLWSRAVI GWSMSPRMTA QLACDALQMA LWRRKRPRNV IVHTDRGGQY CSADYQAQLK RHNLRGMSMA KGCCYDNACV ESFFHSLKVE CIHGEHFISR EIMRATVFNY IECDYNRWRR HSWCGGLSPE QFENKNLA
Source	Mammalian cell
Target Names	insF6
Protein Names	Recommended name: Transposase insF for insertion sequence IS3fA
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.