



Recombinant Human Proto-oncogene c-Fos (FOS)

Product Code	CSB-EP008790HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P01100
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMFSGFNADY EASSSRCSSA SPAGDLSLYY HSPADSFSSM GSPVNAQDFC TDLAVSSANF IPTVTAISTS PDLQWLQPA LVSSVAPSQT RAPHFPGVPA PSAGAYSRAG VVKMTMTGGRA QSIGRRGKVE QLSPEEEEEKR RIRRERNKMA AAKCRNRRRE LTDTLQAETD QLEDEKSALQ TEIANLLKEK EKLEFILAAH RPACKIPDDL GFPEEMSVAS LDLTGGLPEV ATPESSEEAFT LPLLNDPEPK PSVEPVKSIS SMELKTEPFD DFLFPASSRP SGSETARVSP DMDLSGSFYA ADWEPLHSGS LGMGPMATEL EPLCTPVVTC TPSCCTAYTSS FVFYTYPEADS FPSCAAHRK GSSSNPSSD SLSSPTLLAL
Source	E.coli
Target Names	FOS
Protein Names	Recommended name: Proto-oncogene c-Fos Alternative name(s): Cellular oncogene fos G0/G1 switch regulatory protein 7
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. In some cases, expression of the FOS gene has also been associated with apoptotic cell death.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.