



Recombinant *Saccharomyces cerevisiae* 60S ribosomal protein L1-A (RPL1A)

Product Code	CSB-YP317395SVG
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CX43
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SKITSSQVR EHVKELLKYS NETKKRNFL E TVELQVGLKN YDPQRDKRFS GSLKLPNCPR PNMSICIFGD AFDVDRAKSC GVDAMSVDDL KKLNKNKKLI KKLSKKYNAF IASEVLIQV PRLGQPQLSK AGKFPTPVSH NDDLYGKVT VRSTIKFQLK KVLCLAVAVG NVEMEEDVLV NQILMSVNFV VSLKKNWQN VGSLVVKSSM GPAFRLY
Source	Yeast
Target Names	RPL1A
Protein Names	Recommended name: 60S ribosomal protein L1-A Alternative name(s): L10a
Expression Region	2-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.