



Recombinant Escherichia coli Transposase insH for insertion sequence element IS5-12 (JW1882)

Product Code	CSB-EP317611ENV-B
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CE59
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MFVIWSHRTG FIMSHQLTFA DSEFSSKRRQ TRKEIFLSRM EQILPWQNMV EVIIEPFYPKA GNGRRPYPLE TMLRIHCMQH WYNLSDGAME DALYEIASMR LFARLSLDSA LPDRTTIMNF RHLLEQHQLA RQLFKTINRW LAEAGVMMTQ GTLVDATIIE APSSTKNKEQ QRDPEMHQTK KGNQWHFGMK AHIGVDAKSG LTHSLVTTAA NEHDLNQLGN LLHGEEQFVS ADAGYQGAPQ REELAEVDVD WLIAERPGKV RTLKQHPRKN KTAINIEYMK ASIRARVEHP FRIIKRQFGF VKARYKGLLK NDNQLAMLFT LANLFRADQM IRQWERSH
Source	E.coli
Target Names	JW1882
Protein Names	Recommended name: Transposase insH for insertion sequence element IS5-12
Expression Region	1-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.