



# Recombinant Pig Cholecystokinin (CCK)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP004771PI                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                                                           |
| <b>Uniprot No.</b>       | P01356                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Sus scrofa (Pig)                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | QPVPPADSAV PGAQEEEAHR RQLRAVQKVD GESRAHLGAL LARYIQQARK<br>APSGRVSMIK NLQSLDPSHR ISDRDYMGMW DFGRRSAEEY EYTS                                                                                                                                                                                                                                                                                                  |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Target Names</b>      | CCK                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Protein Names</b>     | Recommended name: Cholecystokinin Short name= CCK Cleaved into the following 10 chains: 1. Cholecystokinin-58 Short name= 2. CCK58 3. Cholecystokinin-58 desnonopeptide Alternative name(s): (1-49)-CCK58 Cholecystokinin-39                                                                                                                                                                                |
| <b>Expression Region</b> | 21-114                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Target Details</b>    | Cholecystokinin is a brain/gut peptide. In the gut, it induces the release of pancreatic enzymes and the contraction of the gallbladder. In the brain, its physiologic role is unclear. The cholecystokinin pro-hormone is processed by endo- and exo-proteolytic cleavages.                                                                                                                                |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                             |