



Recombinant Human Tripartite motif-containing protein 49 (TRIM49)

Product Code	CSB-BP317654HU
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0CI25
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNSGILQVFQ GELICPLCMN YFIDPVTIDC GHSFCRPCFY LNWQDIPFLV QCSECTKSTE QINLKTNIHL KKMASLARKV SLWFLSSEE QMCGTHRETK KIFCEVDRSL LCLLCSSSQE HRYHRHRPIE WAAEEHREKL LQKMQLWEK ACENHRNLNV ETTRTRCWKD YVNLRLAIR AEYQKMPAFH HEEKHNLNEM LKKKGKEIFH RLHLSKAKMA HRMEILRGM Y EELNEMCHKP DVELLQAFGD ILHRSESVLL HMPQPLNPEL SAGPITGLRD RLNQFRVHIT LHHEEANNDI FLYEILRSMC IGCDHQDVPY FTATPRSFLA WGVQFTFTSGK YYWEVHVGDS WNWAFGVCNM YRKEKNQNEK IDGKAGLFL GCVKNDIQCS LFTTSPLMLQ YIPKPTSRVG LFLDCEAKTV SFVDVNQSSL IYTIPNCSFS PPLRPIFCCI HF
Source	Baculovirus
Target Names	TRIM49
Protein Names	Recommended name: Tripartite motif-containing protein 49 Alternative name(s): RING finger protein 18 Testis-specific RING-finger protein
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.