



Recombinant *Vaejovis mexicanus smithi* Potassium channel toxin alpha-KTx 21.1

Product Code	CSB-BP317987VAK
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P0DJ31
Product Type	Recombinant Protein
Immunogen Species	<i>Vaejovis mexicanus smithi</i> (Scorpion) (<i>Vaejovis smithi</i>)
Purity	>85% (SDS-PAGE)
Sequence	AAAI SCV GSP ECPPK CRAQG CKNGKCMNRK CKCY YC
Source	Baculovirus
Protein Names	Recommended name: Potassium channel toxin alpha-KTx 21.1 Alternative name(s): Toxin Vm24
Expression Region	1-36
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.