



Recombinant Escherichia coli Transcriptional repressor pifC (pifC)

Product Code	CSB-BP318003ENV
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P10030
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MLSQLNLRFH KKLIEALKTR AGRENTSVNA LAERFLDDGL KTVAPGDGYF QLIADPEATV RQLYRHIILG QTFGTSALSR DELRFVLVHV REAFLRGHNR LATLPALDTL LDITGNLLAW QVEHDRPVDG HYLKGFRLA GKNWTEEFEA FRAALRPVVD QMYAEHLLRP LESDCFGLAE VPDVLAEIF TLPRLKAVFP LMLRGLDWNT EQARTLAQEL RPVISAVTET IEAGTLRLEI RVDGQHPGER PGAWYTTT PRL HLLITGQDFV VPYGWEALSE LLGLFTLYAR HPEALTHGHQ GERVMFSPPG NVTPEGFFGI DGLRIFMPAE AFETLVRELA TRCQEGPLAE ALTGLRCLYG DL
Source	Baculovirus
Target Names	pifC
Protein Names	Recommended name: Transcriptional repressor pifC Short name= Protein C
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.