



# Recombinant Escherichia coli Transcriptional repressor pifC (pifC)

<b>Product Code</b>	CSB-MP318003ENV
<b>Storage</b>	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
<b>Uniprot No.</b>	P10030
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           MLSQLNLRFH KKLIEALKTR AGRENTSVNA LAERFLDDGL KTVAPGDGYF            QLIADPEATV RQLYRHIILG QTFGTSALSR DELRFVLVHV REAFLRGHNR            LATLPALDTL LDITGNLLAW QVEHDRPVDG HYLKGFRLA GKNWTEEFEA            FRAALRPVVD QMYAEHLLRP LESDCFGLAE VPDVLAEIF TLPRLKAVFP            LMLRGLDWNT EQARTLAQEL RPVISAVTET IEAGTLRLEI RVDGQHPGER            PGAWYTTT PRL HLLITGQDFV VPYGWEALSE LLGLFTLYAR HPEALTHGHQ            GERVMFSPPG NVTPEGFFGI DGLRIFMPAE AFETLVRELA TRCQEGPLAE            ALTGLRCLYG DL         </p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	pifC
<b>Protein Names</b>	Recommended name: Transcriptional repressor pifC Short name= Protein C
<b>Expression Region</b>	1-362
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.