



Recombinant Escherichia coli HTH-type transcriptional regulator LeuO (leuO)

Product Code	CSB-MP318016ENV
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P10151
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MPEVQTDHPE TAELSKPQLR MVDLNLLTVF DAVMQEQNIT RAAHVLGMSQ PAVSNVAVARL KVMFNDEL FV RYGRGIQPTA RAFQLFGSVR QALQLVQNEL PGSGFEPASS ERVFHLCVCS PLDSILTSQI YNHIEQIAPN IHVMFKSSLN QNTEHQLRYQ ETEFVISYED FHRPEFTSVP LFKDEMVLVA SKNHPTIKGP LLKHDVYNEQ HAAVSLDRFA SFSQPWYDTV DKQASIAYQG MAMMSVLSV SQTHLVAIAP RWLAEEFAES LELQVLPLPL KQNSRTCYLS WHEAAGRDKG HQMEEQLVS ICKR
Source	Mammalian cell
Target Names	leuO
Protein Names	Recommended name: HTH-type transcriptional regulator LeuO
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.