



Recombinant Rat Protein S100-G (S100g)

Product Code	CSB-MP020644RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P02634
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	SAKKSPEEM KSIFQKYAAK EGDPNQLSKE ELKLLIQSEF PSLLKASSTL DNLFKELDKN GDGEVSYEEF EVFFKLSQ
Source	Mammalian cell
Target Names	S100g
Protein Names	Recommended name: Protein S100-G Alternative name(s): 9 kDa CaBP Calbindin-D9k Cholecalciferol-binding protein S100 calcium-binding protein G Vitamin D- dependent calcium-binding protein, intestinal Short name= CABP
Expression Region	2-79
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes calbindin D9K, a vitamin D-dependent calcium-binding protein. This cytosolic protein belongs to a family of calcium-binding proteins that includes calmodulin, parvalbumin, troponin C, and S100 protein. In the intestine, the protein is vitamin D-dependent and its expression correlates with calcium transport activity. The protein may increase Ca ²⁺ absorption by buffering Ca ²⁺ in the cytoplasm and increase ATP-dependent Ca ²⁺ transport in duodenal basolateral membrane vesicles.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.