



Recombinant *Saccharomyces cerevisiae* Mitochondrial transcription factor 2 (MTF2)

Product Code	CSB-BP318084SVG
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P10849
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	IHCIH RCLLNEANLK DRKTHNVERV SNEKTFEQAL EEERKVFGE FEAGARVENM RHTNASKIID KYYNGLQDNS EGTSVKKEKI VFNHSQRAQR KLPNKDHEFL KETAGNDYVY ERAEPSAIST KTISEQTRTL LEKIFDEDNS INKSNRELLN LNLKRGSGME ALRQVAHSN VKFSEEVMQE IGNKIRYQTT LDQVLEPHID YLREAVKSDY DLLRYLKQSL DIYKRNKDL ELKMNAESSN IFEDIRSACI NKPAELPKPL AMTLPYIIVK SLRLGDFDFP ADRKYTLISY VYNECKNNMD ASLYLTICNV DFYNLLVQLL WENFQEIRYL RRVVTEMSVN GVIGNIETVD ILDKIVKEMR SLNEDVFLEA GEQLSADEEV SSSANKIVNV GVLWNKDTNN DLLIVENYK SLKKNLTRDR
Source	Baculovirus
Target Names	MTF2
Protein Names	Recommended name: Mitochondrial transcription factor 2 Alternative name(s): Protein NAM1
Expression Region	16-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.