



Recombinant *Rhizobium meliloti* Transcriptional regulatory protein fixJ (fixJ)

Product Code	CSB-EP318092RKU
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P10958
Product Type	Recombinant Protein
Immunogen Species	<i>Rhizobium meliloti</i> (strain 1021) (<i>Ensifer meliloti</i>) (<i>Sinorhizobium meliloti</i>)
Purity	>85% (SDS-PAGE)
Sequence	MTDYTVHIVD DEEPVRKSLA FMLTMNGFAV KMHQSAEAFI AFAPDVRNGV LVTDLRMPDM SGVELLRNLG DLKINIPSIV ITGHGDVPM AVEAMKAGAVD FIEKPFEDTV IIEAIERASE HLVAEADVD DANDIRARLQ TLSERERQVL SAVVAGLPNK SIAYDLDISP RTVEVHRANV MAKMKAKSLP HLVRMALAGG FGPS
Source	<i>E.coli</i>
Target Names	fixJ
Protein Names	Recommended name: Transcriptional regulatory protein fixJ
Expression Region	1-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.