



Recombinant Escherichia coli F17 fimbrial protein (F17a-A)

Product Code	CSB-EP318126ENL-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11312
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	YDGKITFNG KVVVDQTC SVT TESKNLTVKL PTVSANS LAS SGKVVGLTPF TILLEGCNTP AVTGAQNVNA YFEPNANTDY TTGNLNTAS SGASNVQIQL LNADGVKAIK LGQAAAAQSV DTVAIN DANV TLR YNAQYYA TGVATAGDVT STVNYTIAYQ
Source	E.coli
Target Names	F17a-A
Protein Names	Recommended name: F17 fimbrial protein Alternative name(s): F17 fimbrial adhesive antigen F17 pilin
Expression Region	22-180
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.