



Recombinant *Drosophila melanogaster* Serine/threonine-protein phosphatase PP-Y (PpY-55A)

Product Code	CSB-MP318156DLU
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P11612
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAVLTTHEID CIKELTSLN GSECTLKEEL IERLIQQTRE VIKWQPMILLE LQAPVNICGD IHGQFTDLLR IFKACGFPPK ANYLFLGDYV DRGKQSLETI CLLFAYKVKY PLNFFLLRGN HESASINKIY GFYDEIKRRH TVRLWHSFTD CFNWLPVAAL VGERIFCCHG GLSPSLRNLQ QINHIQRPTD IPDEGIMCDL LWADLNHTTK GWGHNDRGVS FTFDKVIVRD FLKAFDLQLM VRAHEVVEDG YEFFANRQLV TVFSAPNYCG MMNNAGGVMS VSTDLICSFV IILPCHKYKM IATDANQMPT NEEE
Source	Mammalian cell
Target Names	PpY-55A
Protein Names	Recommended name: Serine/threonine-protein phosphatase PP-Y EC= 3.1.3.16
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.