



Recombinant Escherichia coli 6-phospho-beta-glucosidase BglB (bglB)

Product Code	CSB-MP318190ENV
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P11988
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKAFPETFLW GGATAANQVE GAWQEDGKGI STSDLQPHGV MGKMEPRILG KENIKDVAID FYHRYPEDIA LFAEMGFTCL RISIAWARIF PQGDEVEPNE AGLAFYDRLF DEMAQAGIKP LVTLSHYEMP YGLVKNYGGW ANRAVIDHFE HYARTVFTRY QHKVALWLTF NEINMSLHAP FTGVGLAEES GAEVYQAIH HQLVASARAV KACHSLLPEA KIGNMLLGGL VYPLTCQPQD MLQAMEENRR WMFFGDVQAR GQYPGYMQRFRDHNITIEM TESDAEDLKH TVDFISFSYY MTGCVSHDES INKNAQGNIL NMIPNPHLKS SEWGWQIDPV GLRVLLNTLW DRYQKPLFIV ENGLGAKDSV EADGSIQDDY RIAYLNDHLV QVNEAIADGV DIMGYTSWGP IDLVSASHSQ MSKRYGFIYV DRDDNGEGSL TRTRKKSFGW YAEVIKTRGL SLKKITIKAP
Source	Mammalian cell
Target Names	bglB
Protein Names	Recommended name: 6-phospho-beta-glucosidase BglB EC= 3.2.1.86
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.