



# Recombinant *Xenopus laevis* Serine/threonine-protein kinase mos (mos)

<b>Product Code</b>	CSB-EP318302XBE-B
<b>Storage</b>	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
<b>Uniprot No.</b>	P12965
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus laevis</i> (African clawed frog)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPSPIPVERF LPRDLSPSID LRPCSSPLEL SHRKLPGGLP ACSGRPRLLP PRLAWCSIDW EQVLLLEPLG SGGFGSVYRA TYRGETVALK KVKRSTKNSL ASRQSFVAEL NAARLRHPHV VRVVAASASC PGDPGCPGTI IMEYTGTTGL HQRIYGRSPP LGAEICMRYA RHVADGLRFL HRDGVVHLDL KPANVLLAPG DLCKIGDFGC SQRLREGDEA AGGEPCTQL RHVGGTYTHR APELLKGEPV TAKADIYSFA ITLWQMVSRE LPYTGDRQCV LYAVVAYDLR PEMGPLFSHT EEGRAARTIV QSCWAARPQE RPNAEQLLER LEQECAMCTG GPPSCSPESN APPPLGTGL
<b>Source</b>	<i>E.coli</i>
<b>Target Names</b>	mos
<b>Protein Names</b>	Recommended name: Serine/threonine-protein kinase mos EC= 2.7.11.1 Alternative name(s): pp39-mos
<b>Expression Region</b>	1-359
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.