



# Recombinant Escherichia coli Anaerobic glycerol-3-phosphate dehydrogenase subunit B (glpB)

<b>Product Code</b>	CSB-BP318310ENV
<b>Storage</b>	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
<b>Uniprot No.</b>	P13033
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p> MRFDTVIMGG GLAGLLCGLQ LQKHGLRCAI VTRGQSALHF SSGSLDLLSH  LPDGQPVTDI HSGLESLRQQ APAHPYSLE PQRVLDLACQ AQALIAESGA  QLQGSVELAH QRVTP LGTLR STWLSSPEVP VWPLPAKKIC VVGISGLMDF  QAHLAAASLR ELGLAVETA E IELPELDVLR NNATEFRAVN IARFLDNEEN  WPLLLDALIP VANTCEMILM PACFGLADDK LWRWLNEKLP CSLMLLPTLP  PSVLGIRLQN QLQRQFVRQG GVWMPGDEVK KVTCKNGVVN EIWTRNHADI  PLRPRFAVLA SGSFFSGGLV AERNGIREPI LGLDVLQTAT RGEWYKGDFF  APQPWQQFGV TTDETLRPSQ AGQTIENTLFA IGSVLGGFDP IAQGC GGGVC  AVSALHAAQQ IAQRAGGQQ </p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	glpB
<b>Protein Names</b>	Recommended name: Anaerobic glycerol-3-phosphate dehydrogenase subunit B Short name= Anaerobic G-3-P dehydrogenase subunit B Short name= Anaerobic G3Pdhase B EC= 1.1.5.3
<b>Expression Region</b>	1-419
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.