



Recombinant Rhizobium meliloti Nitrogen fixation regulation protein fixK (fixK)

Product Code	CSB-BP318344RKU
Storage	The shelf life of liquid-form form is 6 months around at -20°C/-80°C. The shelf life of lyophilized form is 12 months around at -20°C/-80°C.
Uniprot No.	P13295
Product Type	Recombinant Protein
Immunogen Species	Rhizobium meliloti (strain 1021) (Ensifer meliloti) (Sinorhizobium meliloti)
Purity	>85% (SDS-PAGE)
Sequence	MYAAAQAKPQ SIEVEHLGPA PMSGPRLVAT YKPGREIYAQ GDLNDKCYQV STGAVRIYRL LSDGRRQVVS FHLPGEMFGF EAGSNHSFFA EAITETTLAI FGRRNMQERS RELLALALTG MARAQHLLV IGRQCAVERI AAFLVDLCER QGGGRQLRLP MSRQDIADYL GLTIETVSRV VTKLKERSLI ALRDARTIDI MKPEALRSLC N
Source	Baculovirus
Target Names	fixK
Protein Names	Recommended name: Nitrogen fixation regulation protein fixK
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.